**R<sub>x</sub> ONLY**

**STERILE R**

## Indications For Use

The CBCII is a closed blood recovery system used post-operatively to collect, filter, and allow for reinfusion of autologous blood.

## Contraindications

- Abnormal renal and/or hepatic function
- Malignant lesions
- Contamination/sepsis
- Fluids unsuitable for reinfusion
- Amniotic fluid
- Bile
- Hemostatic agents
- Excessive hemolysis
- Coagulation disorders
- Potential for air embolism, microembolism, fat embolism

## Description

The sterile CBCII ConstaVac Blood Conservation system is a single patient use disposable kit. Each kit typically contains the following components:

- 1 CBCII system
- A variety of wound drain types, quick disconnect kits and/or other accessories

**NOTE:** For a complete list of accessory information, contact your 3T Medical sales representative. Outside the US, contact your nearest 3T Medical subsidiary.

## User/Patient Safety

### WARNINGS:

- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use.
- Only trained and experienced healthcare professionals should use this equipment.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. 3T Medical, as a manufacturer, does not recommend surgical procedure or technique.
- The CBCII system is used specifically for post-operative reinfusion of autologous blood, not allogenic blood. If complications arise, use appropriate clinical management to monitor coagulation factors related to the reinfusion of autologous blood.
- Upon initial receipt and before use, inspect each component for damage. DO NOT use any component if damage is apparent.
- Use only 3T Medical-approved system components and accessories, unless otherwise specified. Using other electronic accessories may result in increased electromagnetic emissions or decreased electromagnetic immunity of the system.
- DO NOT modify any system component or accessory.
- DO NOT service this equipment.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like this system. Install and place this system into service according to the EMC information contained in this manual. Portable and mobile RF communications equipment can affect the function of this system.

### WARNINGS:

- DO NOT use this equipment in areas in which flammable anesthetics or flammable agents are mixed with air, oxygen or nitrous oxide.
- Upon initial receipt and before use, inspect the package for damage and confirm the integrity of the sterile barrier. DO NOT use any equipment if damage is apparent or the sterile barrier has been compromised.
- DO NOT use the equipment after the expiration date printed on the package. The equipment may not be safe or effective after the expiration date.
- Single patient use only.
- DO NOT reuse, reprocess, or repackage this device. This device is intended for a single use only. This device may not withstand chemical, chemical vapor, or high temperature sterilization reprocessing. Design features may make cleaning difficult. Reuse may create a serious risk of contamination and may compromise the structural integrity of the device resulting in operational failure. Critical product information may be lost if the device is repackaged. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.
- Dispose of properly.
- ALWAYS follow the current local regulations governing the prevention of transmission of bloodborne infections when handling any blood-related product.
- ALWAYS keep the system upright during use. DO NOT turn the system sideways or upside down. Failure to comply may cause the air filter to clog and prevent transfer of blood from the system.

## Introduction

This *Instructions For Use* manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, surgical technologists, and biomedical equipment technicians. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

If additional information, especially safety information, or in-service training is required, contact your 3T Medical sales representative or call 3T Medical customer service. Outside the US, contact your nearest 3T Medical subsidiary.

Trademarks not the property of 3T Medical Systems, Inc. are the property of their respective owners.


ENGLISH (EN)	PORTUGUÊS (PT)
ESPAÑOL (ES)	NORSK (NO)
DEUTSCH (DE)	POLSKI (PL)
FRANÇAIS (FR)	ΕΛΛΗΝΙΚΑ (EL)
ITALIANO (IT)	日本語 (JA)
NEDERLANDS (NL)	中文 (ZH)
SVENSKA (SV)	한국어 (KO)
DANSK (DA)	
SUOMI (FI)	

### Legend


Figure	Caption
<b>A</b>	Features 1 - Y-connector 2 - Luer Connector 3 - Cap 4 - Reservoir Label 5 - Drain Lever 6 - Vacuum Control Dial 7 - Vacuum Indicator 8 - Mounting Clamp 9 - Security Strap 10 - Collection Reservoir (reservoir) 11 - Evacuator Tube 12 - Reservoir/Blood Bag Tube 13 - Slide Clamp 14 - Blood Bag (reinfusion bag)
<b>B</b>	To Prepare/Install the System
<b>C</b>	To Collect Blood Using Vacuum
<b>D</b>	To Transfer/Reinfuse Collected Blood
<b>E</b>	To Collect/Dispose of Blood - No Reinfusion (optional)
<b>F</b>	To Remove/Replace the System (optional)

### Symbol Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

			
Off	Low Vacuum	Medium Vacuum	High Vacuum
			
Refer to instruction manual/booklet	General warning sign	Consult instructions for use	DO NOT reuse
			
Direct Current (DC)	Sterilized by irradiation	Fully Push Down To Drain See Operating Instructions	Keep Upright In Use

**A**

## Instructions

### To Prepare/Install the System - See Figure B

**WARNING:** Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use. See the *User/Patient Safety* section.

**CAUTION:** ALWAYS operate the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.

#### NOTES:

- The entire CBCII device is an applied part as defined by the manufacturer according to the standards listed under *Product Safety Certification* in the *Specifications* section.
- The packaging will include the blood conservation system and optional accessories, including wound drains.
- Keeping the blue cover wrapped around the system is recommended to prevent blood transfer from the wound site to the system.
- Remove and discard the blue sterile cover properly after use.

#### WARNINGS:

- Remove the protective cover from the trocar with care.
- ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of sharps and single-use medical devices.

1. Using a trocar or needle, insert and place the wound drain into the wound site as directed by local hospital protocol.
2. Remove the evacuator tube with the Y-connector from within the sterile blue cover. Trim the Y-connector at the proper transitional section.

**NOTE:** Each transitional section of the Y-connector corresponds to a specific diameter of wound drain tubing.

**WARNING:** Make sure all connections are airtight. DO NOT reinfuse collected blood if bubbles are observed in the reservoir or the evacuator tube. If an air leak occurs, reinfusion of blood is contraindicated.

**CAUTION:** If bone cement is in a hardened state, irrigate the wound to remove any excess cement material. Failure to comply may cause the system to clog.

3. Install the wound drain tubing into the Y-connector. Make sure each connection is sterile and air tight.

#### NOTES:


- Applying surgical tape to the connection between the Y-connector and wound drain tubing is recommended.
- If a tourniquet is used, wait 15 minutes after the tourniquet is released before turning ON the blood conservation system. Waiting will help to eliminate the occlusion of the drain sites by adjacent tissue that might prevent the flow of blood to the system.
- Make sure the system is oriented in an upright position at all times.

**WARNING:** ALWAYS keep the system upright during use. DO NOT turn the system sideways or upside down. Failure to comply may cause the air filter to clog and prevent transfer of blood from the system.

4. Orient the system in an upright position near the patient.

## B

## Instructions

### To Collect Blood Using Vacuum - See Figure C

**WARNING:** Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use. See the *User/Patient Safety* section.

#### CAUTIONS:

- When relocating the patient with the system, make sure the system is secure and protected from physical damage. The reservoir should be oriented toward the patient.
  - Empty the reservoir before moving the patient. Failure to comply may cause the air filter to clog.
1. To relocate the patient with the system, mount the system to the foot board or side rail of the bed. Use the security strap and adjustable mounting clamp to secure the system to the bed before transport.
  2. Turn the vacuum control dial to the I setting (LOW). See *Specifications* for vacuum level values. If the system fails to operate, see the *Troubleshooting* section.
  3. Observe the vacuum indicator is inverted. If the indicator does not invert, see the *Troubleshooting* section.
  4. Record the patient name, start time and the volume level on the reservoir label. Allow blood to collect at the LOW vacuum setting for at least 10 minutes.
  5. After 10 minutes, turn the vacuum control dial to a setting prescribed by the physician.

**NOTE:** To maintain an accurate vacuum level setting, the system should be located at the approximate level of the wound site. The vacuum level is reduced by 0.5 psi [25.85 mmHg] for every 12 inches [30 cm] the system is elevated above the wound site.


6. Monitor the volume of the reservoir and the rate of blood collection directly from the reservoir. If blood collection is minimal, or is not occurring, or blood clots are forming inside the reservoir, see the *Troubleshooting* section.

**WARNING:** ALWAYS monitor the volume of the reservoir and the rate of blood collection. Failure to comply may allow excessive blood loss to occur unnoticed.

**NOTE:** Leave the blood bag banded until ready to transfer collected blood.

7. To stop collecting blood and end the operation of the system, turn the vacuum control dial to the O setting (OFF).

## Instructions

### To Transfer/Reinfuse Collected Blood - See Figure D

**WARNING:** Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use. See the *User/Patient Safety* section.

#### NOTES:

- Transfer and reinfusion of collected blood should be performed after collecting 400 mL of blood or after four hours of blood collection, whichever occurs first.
  - Placing the blood bag approximately 12 inches below the reservoir is recommended to improve blood transfer into the blood bag.
1. Unband the blood bag and tubing. Place the blood bag approximately 12 inches below the reservoir.
  2. Hold the blood bag tubing to form a half loop.

**WARNING:** DO NOT lift the lever. Failure to comply may damage the product and prevent the successful transfer or reinfusion of blood to the patient.

3. Fully push down and hold the lever.
4. When the reservoir is full and blood transfer begins, allow the tubing to hang freely and straighten.

#### NOTES:

- Although the system does not require a minimum amount of blood for reinfusion, less than 70 mL of collected blood may not transfer to the blood bag.
  - If less than 100 mL of blood is collected within four hours, discard the system or use the system for drainage only.
  - Pumping the lever may assist with initiating blood transfer.
  - While holding down the lever, kneading the tube between the reservoir and bag may force the remaining blood into the bag.
5. When blood transfer is complete and the reservoir is empty, release the lever. If blood fails to transfer, see the *Troubleshooting* section.

**NOTE:** The remaining fluid in the reservoir may contain fat that should not be reinfused.

6. Using the slide clamp, close the tubing as close to the blood bag as possible to eliminate the canister vacuum from competing with the gravity drain for the reinfusion.

#### WARNINGS:

- Read and understand the instructions for use supplied with the blood administration kit and filter before use. ALWAYS make sure the initial venting and priming of the kit is performed properly.
- DO NOT reinfuse topical hemostatic agents, wound irrigants, or antibiotics that are not licensed for parenteral use.
- DO NOT reinfuse the collected blood if methylmethacrylate is present in liquid form.
- If necessary, elevate the blood bag to increase the rate of reinfusion. DO NOT use a mechanical device to increase the rate of reinfusion. Failure to comply may increase the risk of air embolism.
- DO NOT mix collected blood with any medication.
- DO NOT add a medication into a flowing blood path.
- If the blood bag is used to store collected blood, make sure the bag is labeled properly to prevent unintentional reinfusion.

7. Reinfuse the collected blood. If collected blood fails to reinfuse from the blood bag to the administration kit, see the *Troubleshooting* section.

#### NOTES:


- Reinfuse the collected blood using a 20 or 40 micron microaggregate filter and blood administration kit.
  - Reinfuse as directed by standard hospital protocol.
8. Repeat the process for additional reinfusions or discontinue using the blood conservation system as directed by standard hospital protocol.

**NOTE:** The system remains a closed system after each reinfusion if the instructions for use are followed precisely.

#### WARNINGS:

- Remove the wound drain from the wound site with care.
  - ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of sharps and single-use medical devices.
9. Remove the blood conservation system and accessories from the patient. Properly dispose of all sharps and biohazard waste. See the *Disposal/Recycle* section.

## Instructions

### To Collect/Dispose of Blood - No Reinfusion (optional)

**WARNING:** Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use. See the *User/Patient Safety* section.

#### NOTES:

- The system may be used to collect, transfer, and dispose of blood in a blood bag as a one time occurrence. However, the system may still be used as a wound drainage receptacle. No reinfusion is possible.
- Alternatively, the blood bag may be removed and discarded before use. The system may then be used to collect blood only. No reinfusion is possible.
- ALWAYS measure and monitor the collected blood directly from the reservoir.

### BLOOD DISPOSAL WITH BLOOD BAG

#### Initial System Use

1. Prepare and install the system. See *To Prepare/Install the System*.
2. Collect blood in the system. See *To Collect Blood Using Vacuum*.
3. Transfer the collected blood into the blood bag. See *To Transfer/Reinfuse Collected Blood*.

#### Blood Bag Removal - See Figure E

1. Use the slide clamp to close the tubing as close to the reservoir as possible.
2. Cut the blood bag tubing about two inches [five cm] away from the reservoir. Use the cap to close the tube.
3. Discard and properly dispose of the blood-filled bag and tubing as directed by standard hospital protocol.

#### Final System Use

1. Collect additional blood in the system as required. See *To Collect Blood Using Vacuum*.
2. To empty the reservoir, place a separate blood collection container under the reservoir tube, remove the cap, and use the slide clamp to open, and then close the tubing. Replace the cap.
3. Discard and properly dispose of the collected blood as directed by standard hospital protocol.

#### WARNINGS:

- Remove the wound drain from the wound site with care.
  - ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of sharps and single-use medical devices.
4. Discontinue using the system as directed by standard hospital protocol. Remove the wound drain(s), tubing, and system from the patient. Properly dispose of sharps and all biohazard waste. See the *Disposal/Recycle* section.

### BLOOD DISPOSAL WITHOUT BLOOD BAG

#### Blood Bag Removal - See Figure E

1. Using the slide clamp, close the tubing as close to the reservoir as possible.
2. Cut the blood bag tubing about two inches [five cm] away from the reservoir. Using the cap, close the tube.
3. Discard and properly dispose of the empty blood bag and tubing.

#### System Use


1. Prepare and install the system. See *To Prepare/Install the System*.
2. Collect blood in the system. See *To Collect Blood Using Vacuum*.
3. To empty the reservoir, place a blood collection receptacle under the reservoir tube, remove the cap, and use the slide clamp to open and close the tubing. Replace the cap as required.
4. Discard and properly dispose of the collected blood as directed by standard hospital protocol.

#### WARNINGS:

- Remove the wound drain from the wound site with care.
  - ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of sharps and single-use medical devices.
5. Discontinue using the system as directed by standard hospital protocol. Remove the wound drain(s), tubing, and system from the patient. Properly dispose of sharps and all biohazard waste. See the *Disposal/Recycle* section.

# E

## Instructions

### To Remove/Replace the System (optional) - See Figure F

**WARNING:** Before using any system component, or any component compatible with this system, read and understand the instructions. Pay special attention to WARNING information. Become familiar with the system components prior to use. See the *User/Patient Safety* section.

1. Use the slide clamps to close the evacuator tube.
2. Disconnect and separate the Luer connector.
3. Remove and properly dispose of the old system as biohazard waste.
4. Obtain a new system. Disconnect the evacuator tube with the Y-connector from the new system and dispose of the tube properly.
5. Connect the new system to the evacuator tube attached to the patient using the Luer connector.
6. Turn the vacuum control dial to a setting prescribed by the physician.
7. Use the slide clamp(s) to open the evacuator tube.

## Disposal/Recycle

### WARNINGS:

- ALWAYS follow the current local regulations governing biohazard waste to safely handle and dispose of sharps and single-use medical devices.
- Handling biohazard waste is potentially dangerous. ALWAYS follow current local regulations governing biohazard waste to safely handle and dispose of surgical waste.
- ALWAYS follow the current local recommendations and/or regulations governing environmental protection and risks associated with recycling or disposing of the equipment at the end of its useful life.
- ALWAYS decontaminate a medical device exposed to infectious material before sending it to a waste treatment facility. If decontamination of the device is not possible, DO NOT remove the batteries from the device.


To comply with Batteries Directive 2006/66/EC, this device has been designed for safe removal of the batteries at the end-of-life by a waste treatment facility. If decontamination is not possible, Batteries Directive 2006/66/EC and Member State regulations allow the disposal of small amounts of portable batteries to landfill and incineration.

To comply with European Community Waste Electrical and Electronic Equipment (WEEE) Directive 2012/19/EU, ALWAYS collect this product separately for recycling. DO NOT dispose of this product as unsorted municipal waste. Contact your local distributor for disposal information.

## F

## Storage and Handling

**CAUTION:** ALWAYS store the equipment within the specified environmental condition values throughout its useful life. See the *Specifications* section.

To ensure the longevity, performance and safety of this equipment, use of the original packaging material is recommended when storing or transporting this equipment.

## Troubleshooting

**WARNING:** DO NOT service this equipment, unless otherwise specified.

**NOTE:** If you require service, contact your 3T Medical sales representative or call 3T Medical customer service. Outside the US, contact your nearest 3T Medical subsidiary.

PROBLEM	CAUSE	SOLUTION
The system fails to work.	A vacuum lock may exist in the reservoir.	Press and hold the lever or disconnect the luer connector to release the vacuum.
	The system is damaged.	If necessary, replace the system. See <i>To Remove/Replace the System</i> .
The vacuum indicator fails to invert.	The vacuum control dial may not be turned on.	Turn the vacuum control dial to the prescribed setting.
	The system may have an air leak.	Inspect all connections. Make all connections air tight. If air leak is suspected, contact a physician. Reinfusion is contraindicated.
	The system air filter may be saturated or the system is damaged.	If necessary, replace the system. See <i>To Remove/Replace the System</i> .
Minimal or no blood is collected by the system.	Wound drain tubing or evacuator tubing may be blocked.	Knead the tubing to stimulate flow. Turn the vacuum control dial to a higher vacuum setting.
	Clots may be in the Y-connector.	Knead the clots through the Y-connector.
	The reservoir is located above the wound site.	Relocate the system at or below the wound site. Turn the vacuum control dial to a higher vacuum setting.
	The system may have an air leak.	Inspect all connections. Make all connections air tight. If air leak is suspected, contact a physician. Reinfusion is contraindicated.
	The wound drain may not be placed correctly.	Contact physician to relocate wound drain.
	The system pre-filter may be saturated or the system is damaged.	If necessary, replace the system. See <i>To Remove/Replace the System</i> .
	Blood clots appear in the reservoir.	A tourniquet was applied at the wound site and the vacuum control dial was turned on too soon after the tourniquet was removed.
Blood fails to transfer from the reservoir to the blood bag.	A vacuum lock may exist in the reservoir.	Pump the drain lever or hold down drain lever and knead the tubing to stimulate flow.
	The blood bag is located above the reservoir.	Relocate the blood bag approximately 12 inches below the reservoir.
	The system air filter may be saturated or the system is damaged.	 <p><b>WARNING:</b> Inserting a needle into the vacuum indicator will allow unfiltered air into the reservoir. Reinfusion of the blood is contraindicated.</p> <p>Insert a needle into the vacuum indicator. Press the drain lever to transfer the blood from the reservoir. Replace the system. See <i>To Remove/Replace the System</i>.</p>

PROBLEM	CAUSE	SOLUTION
Blood fails to transfer from the blood bag to the administration kit.	The administration kit filter may not be primed or is air locked.	Manually squeeze the blood bag until the blood begins to flow to the kit.
Sporadic electrical interference is experienced.	Electrical noise is present.	Turn off all electrical equipment not in use in the operating room.
		Relocate electrical equipment; increase spatial distance.
		Plug operating room equipment into different operating room outlets.

## Specifications

<b>Model:</b>	CBCII ConstaVac Blood Conservation System	
<b>REF:</b>	0225-XXX-XXX Series	
<b>European Conformity:</b>	 <b>0459</b>	REF 0225-028-000E and REF 0225-028-626 only
<b>Dimensions:</b>	Height: 229 mm [9.0 inch] Width: 125 mm [4.9 inch] Depth: 165 mm [6.5 inch]	
<b>Mass:</b>	1.1 kg [2.5 lb]	
<b>Volume:</b>	Reservoir: 800 mL + 10 mL Blood Bag: 800 mL + 30 mL	
<b>Evacuator Tube:</b>	6 feet [183 cm] length	
<b>Filters:</b>	Short Pre-filter: 200 micron Air filter: 0.45 micron	
<b>Vacuum Level Settings:</b>	Low: I 25 mmHg [0.5 psi] Medium: II 50 mmHg [1.0 psi] High: III < 100 mmHg [< 1.9 psi]	
<b>Power Supply:</b>	Internally powered 3.0 V  , two AA alkaline batteries factory installed Duration: > 24 hours at high vacuum (level III on the vacuum control dial)	
<b>Ingress Protection:</b>	IPX0	
<b>Equipment Type:</b>		Type BF Applied Part

**NOTE:** The entire CBCII device is an applied part as defined by the manufacturer according to the standards listed under *Product Safety Certification* in the *Specifications* section.

## Specifications (continued)

**Product Safety Certification:**

**Canadian Standards Association (CSA) International**

C US

**Canadian Standards Association (CSA)**

CAN/CSA-C22.2 No. 601.1-M90, *Medical Electrical Equipment – Part 1: General Requirements for Safety*

**Underwriters Laboratories (UL)**

UL 60601-1, *Medical Electrical Equipment, Part 1: General Requirements for Safety – First Edition*; Revisions through and including April 26, 2006

For REF 0225-028-000E and REF 0225-028-626 only:

**International Electrotechnical Commission (IEC)**

IEC 60601-1:2005, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* Corrigendum 1 (2006); Corrigendum 2 (2007)

IEC 60601-1:1988, *Medical Electrical Equipment – Part 1: General Requirements for Safety - Second Edition*; Amendment 1 (1991); Amendment 2 (1995); Corrigendum 1 (1995)

**Canadian Standards Association (CSA)**

CAN/CSA-C22.2 No. 60601-1:08, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*

**American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)**

ANSI/AAMI ES60601-1:2005, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* Consolidated Reprint (2009); Amendment 2 (2010)


**European Committee for Electrotechnical Standardization (CENELEC)**

EN 60601-1:2006, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011)


**Environmental Conditions:**

Temperature Limitation:


**Operation**



**Storage and Transportation**

Humidity Limitation:

Atmospheric Pressure Limitation:

## Specifications

**Guidance and Manufacturer's Declaration - Electromagnetic Emissions**

The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series is intended for use in the electromagnetic environment specified below. The customer or the user of the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	Not Applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not Applicable	Not Applicable

**Guidance and Manufacturer's Declaration - Electromagnetic Immunity**

The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series is intended for use in the electromagnetic environment specified below. The customer or the user of the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3.5/V_t] \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  Where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m) Interference may occur in the vicinity of equipment marked with the following symbol:  (Non-ionizing electromagnetic radiation)

NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Specifications (continued)

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series is intended for use in the electromagnetic environment specified below. The customer or the user of the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD)  IEC 61000-4-2	±6 kV contact  ±8 kV air	±2 kV, ±4 kV, ±6 kV contact  ±2 kV, ±4 kV, ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst  IEC 61000-4-4	±2 kV for power supply lines  ±1 kV for input/output lines	Not Applicable	Not Applicable
Surge  IEC 61000-4-5	±1 kV differential mode  ±2 kV common mode	Not Applicable	Not Applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines  IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 seconds	Not Applicable	Not Applicable
NOTE: $U_T$ is the alternating current mains voltage prior to application of the test level.			

## Specifications

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series			
The CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CBCII ConstaVac Blood Conservation System REF 0225-XXX-XXX Series as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter  W	Separation distance according to frequency of transmitter  m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = [3.5 \sqrt{V_T}] \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	Not Applicable	0.12	0.23
0.1	Not Applicable	0.38	0.73
1	Not Applicable	1.2	2.3
10	Not Applicable	3.8	7.3
100	Not Applicable	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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C US

### Canadian Standards Association (CSA)

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### Underwriters Laboratories (UL)

UL 60601-1, *Medical Electrical Equipment, Part 1: General Requirements for Safety – First Edition; Revisions through and including April 26, 2006*

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IEC 60601-1:1988, *Medical Electrical Equipment – Part 1: General Requirements for Safety – Second Edition; Amendment 1 (1991); Amendment 2 (1995); Corrigendum 1 (1995)*

### Canadian Standards Association (CSA)

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### American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)

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### European Committee for Electrotechnical Standardization (CENELEC)

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Canadian Standards Association (CSA) International

C US

**Canadian Standards Association (CSA)**  
CAN/CSA-C22.2 No. 601.1-M90, *Medical Electrical Equipment – Part 1: General Requirements for Safety*

**Underwriters Laboratories (UL)**  
UL 60601-1, *Medical Electrical Equipment, Part 1: General Requirements for Safety – First Edition*; Revisions through and including April 26, 2006

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**International Electrotechnical Commission (IEC)**  
IEC 60601-1:2005, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* Corrigendum 1 (2006); Corrigendum 2 (2007)

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**Canadian Standards Association (CSA)**  
CAN/CSA-C22.2 No. 60601-1:08, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance*

**American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)**  
ANSI/AAMI ES60601-1:2005, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* Consolidated Reprint (2009); Amendment 2 (2010)

**European Committee for Electrotechnical Standardization (CENELEC)**  
EN 60601-1:2006, *Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance* IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011)

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Canadian Standards Association (CSA)  
CAN/CSA-C22.2 No. 601.1-M90, Medical Electrical Equipment — Part 1: General Requirements for Safety

Underwriters Laboratories (UL)  
UL 60601-1, Medical Electrical Equipment, Part 1: General Requirements for Safety — First Edition; Revisions through and including April 26, 2006

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International Electrotechnical Commission (IEC)  
IEC 60601-1:2005, Medical Electrical Equipment — Part 1: General Requirements for Basic Safety and Essential Performance Corrigendum 1 (2006); Corrigendum 2 (2007)

IEC 60601-1:1988, Medical Electrical Equipment — Part 1: General Requirements for Safety - Second Edition; Amendment 1 (1991); Amendment 2 (1995); Corrigendum 1 (1995)

Canadian Standards Association (CSA)  
CAN/CSA-C22.2 No. 60601-1:08, Medical Electrical Equipment — Part 1: General Requirements for Basic Safety and Essential Performance

American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)  
ANSI/AAMI ES60601-1:2005, Medical Electrical Equipment — Part 1: General Requirements for Basic Safety and Essential Performance Consolidated Reprint (2009); Amendment 2 (2010)



European Committee for Electrotechnical Standardization (CENELEC)  
EN 60601-1:2006, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011)

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


C US

Canadian Standards Association (CSA)  
CAN/CSA-C22.2 No. 601.1-M90, Medical Electrical Equipment —  
Part 1: General Requirements for Safety

Underwriters Laboratories (UL)  
UL 60601-1, Medical Electrical Equipment, Part 1: General  
Requirements for Safety — First Edition; Revisions through and  
including April 26, 2006

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International Electrotechnical Commission (IEC)  
IEC 60601-1:2005, Medical Electrical Equipment — Part 1: General  
Requirements for Basic Safety and Essential Performance Corrigendum  
1 (2006); Corrigendum 2 (2007)

IEC 60601-1:1988, Medical Electrical Equipment — Part 1: General  
Requirements for Safety - Second Edition; Amendment 1 (1991);  
Amendment 2 (1995); Corrigendum 1 (1995)

Canadian Standards Association (CSA)  
CAN/CSA-C22.2 No. 60601-1:08, Medical Electrical Equipment  
— Part 1: General Requirements for Basic Safety and Essential  
Performance

American National Standards Institute (ANSI)/Association for the  
Advancement of Medical Instrumentation (AAMI)  
ANSI/AAMI ES60601-1:2005, Medical Electrical Equipment — Part  
1: General Requirements for Basic Safety and Essential Performance  
Consolidated Reprint (2009); Amendment 2 (2010)

European Committee for Electrotechnical Standardization (CENELEC)  
EN 60601-1:2006, Medical Electrical Equipment — Part 1: General  
Requirements for Basic Safety and Essential Performance IEC  
Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC  
Corrigendum (2010); CENELEC Amendment A11 (2011)

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

Canadian Standards Association (CSA)  
 CAN/CSA-C22.2 No. 601.1-M90, Medical Electrical Equipment – Part 1: General Requirements for Safety Underwriters Laboratories (UL)  
 UL 60601-1, Medical Electrical Equipment, Part 1: General Requirements for Safety – First Edition; Revisions through and including April 26, 2006

REF 0225-028-000E  
 및 REF 0225-028-626  
 전용:

International Electrotechnical Commission (IEC)  
 IEC 60601-1:2005, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance Corrigendum 1 (2006); Corrigendum 2 (2007)

IEC 60601-1:1988, Medical Electrical Equipment – Part 1: General Requirements for Safety – Second Edition; Amendment 1 (1991); Amendment 2 (1995); Corrigendum 1 (1995)

Canadian Standards Association (CSA)  
 CAN/CSA-C22.2 No. 60601-1:08, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance

American National Standards Institute (ANSI)/Association for the Advancement of Medical Instrumentation (AAMI)  
 ANSI/AAMI ES60601-1:2005, Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance Consolidated Reprint (2009); Amendment 2 (2010)



European Committee for Electrotechnical Standardization (CENELEC)  
 EN 60601-1:2006, Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance IEC Corrigendum 1 (2006); IEC Corrigendum 2 (2007); CENELEC Corrigendum (2010); CENELEC Amendment A11 (2011)

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