
Colder® Insulated Connector Hose

Operations Manual
REF 8001-064-135

EN NL
CS NO
DA PL
DE PT
ES RU
FR SV
IT TR
KO ZH

8001-009-202 Rev A

2025-09

EN

Symbols

	Refer to instruction manual/booklet
	Consult instructions for use
	General warning
	Caution
	Catalogue number
	Batch code
	Unique device identifier
	Manufacturer
	Date of manufacture
	For US Patents see www.3tmedical.com/patents
	CAUTION - US Federal Law restricts this device to sale by or on the order of a licensed physician.
	In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.
	Temperature controller

Introduction

This manual assists you with the operation or maintenance of your 3T Medical product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.

Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- 3T Medical continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact 3T Medical Customer Service or Technical Support at (734)525-9580.

Product description

The **Colder®** Insulated Connector hose provides an interface from the thermal transfer devices (wraps or blankets) with **Colder®** Connector connections to the **Altrix** Precision Temperature Management System (referred to as the controller).

Expected life

The **Colder®** Insulated Connector Hose has a one year use expected life under normal use conditions.

Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

Specifications

Patient facing material	EPDM (Ethylene Propylene DI-Monomer Sponge Extrusion Flame Retardant Grey Flexible PVC Adhesive Lined, Polyolefin Polypropylene)	
Length	10 ft	3.05 m
Compatible controller	Model 8001 Altrix Precision Temperature Management System	See Model 8001 Operations Manual for cleaning, disinfection, and more information.

3T Medical reserves the right to change specifications without notice.

European REACH

In accordance with the European REACH regulation and other environmental regulatory requirements the components that contain declarable substances are listed.

Description	Number	Substance of very high concern (SVHC) chemical name
Colder® Insulated Connector Hose	8001-064-135	bis (2-ethylhexyl) phthalate

Contact information

Contact 3T Medical Customer Service or Technical Support at: (734)525-9580

3T Medical Systems, Inc.
236 Fernwood Avenue
Edison, NJ 08837
USA

Note - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see <https://3tmedical.com/technical-docs>

Have the lot batch code (A) of your 3T Medical product available when calling 3T Medical Customer Service or Technical Support. Include the serial number in all written communication.

Setup

- Connect the hose (A) to the controller hose connection port (Figure 1).


Figure 1 –

- Connect the hose end (B) to the thermal transfer device (C).

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